


The Gallery Menu ("View" Tab)



1. Overview & Purpose

The **Gallery Menu** (accessed via the "View" tab) acts as the central review hub for the clinical session. It is the bridge between **Data Capture** (Camera) and **Data Utilization** (Reporting & Cloud Sync). Its primary purpose is to allow the practitioner to verify the quality and completeness of the exam data before the patient leaves.

Unlike a standard file viewer, the Gallery organizes media **anatomically** and **contextually**, grouping images into pairs that mirror the clinical workflow. It also serves as the secure gateway for generating reports and syncing data to the cloud.

2. User Interface Structure

The UI follows a strict hierarchy designed for clinical efficiency:

- **Top Navigation (Tabs):**
 - **Upper Arch / Lower Arch:** These tabs display the guided scan sequence. Images are strictly ordered by sequence number (e.g., Sequence #1, #2, #3).
 - **Other:** Holds non-sequence media like full-face photos or video recordings.
 - **The "Paired Sequence" Cards:**
 - Images are not shown as a loose list. They are grouped into **Cards**.
 - Each Card represents a specific anatomical site and contains two slots:
 1. **RGB (Left):** The standard visual light image.
 2. **Fluorescence (Right):** The corresponding QLF diagnostic image.
 - **Logic:** The system automatically pairs images based on their SequenceNumber. If an image is missing (e.g. RGB mode capture hasn't been done yet and only fluorescence mode is captured), a "Waiting for Capture" placeholder is shown, acting as a visual checklist for the doctor.
 - **Action Bar (Footer):**
 - **Discard Pair (Red):** Allows immediate deletion of blurry or incorrect images (deletes from both disk and database).
 - **Get Report (Green):** Triggers the analysis engine to generate a diagnostic PDF.
 - **Save and Sync (Purple):** Encrypts and uploads all pending data to the cloud.
-

3. Data & Security Architecture

The Gallery is the visual interface for the application's **Zero-Trust Storage System**.

- **Secure Storage:** All images displayed in the Gallery are stored on the local disk as **AES-256 Encrypted** files. They cannot be opened by standard Windows tools.
 - **Decrypt-on-Demand:** When the Gallery loads, it reads metadata from the secure SQLite database (AppDbContext). It then decrypts the image file **directly into memory (RAM)** to display the thumbnail. The decrypted version is never written back to the hard drive, ensuring maximum security.
 - **Persistence:** The Gallery's state is persistent. If the app is closed and reopened, the Gallery reloads the patient's exact state from the secure database, ensuring no data loss.
-

4. Workflows & Features

A. Real-Time Updates:

The Gallery listens for global application events. When a new image is captured in the Camera tab, it instantly appears in the Gallery without needing to refresh, providing immediate feedback.

B. Report Generation:

Clicking "Get Report" initiates the analysis pipeline:

1. It aggregates all Fluorescence images from the session.
2. It runs computer vision analysis to score caries and plaque.
3. It generates a PDF report, **encrypts it**, and saves it securely as a new "Report" item in the patient's media record.

C. Cloud Synchronization:

Clicking "Save and Sync" triggers the SyncManager:

1. **Patient Data:** Syncs the patient's demographic record to the cloud.
2. **Media Upload:** Checks the local database for any items (Images, Videos, Reports) where S3Url is null.
3. **Secure Transfer:** It decrypts these files in memory and uploads them to AWS S3 using a secure, patient-centric folder structure: {GlobalPatientId}/{ClinicID}/{FileName}.
4. **Completion:** Updates the local database with the new Cloud URLs to mark the sync as complete.